
Labor

N022.18, N026.18

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
Labor

N022.18, N026.18

**TRANSLATION OF ORIGINAL
OPERATING INSTRUCTION
ENGLISH**

EN

DIAPHRAGM PUMP

Notice!

Before operating the pump and accessories, read and observe the operating and installation instructions as well as the safety information!

Contents

1	About this document	3
1.1	Using the operating instructions.....	3
1.2	Exclusion of liability.....	3
1.3	Symbols and markings.....	4
2	Safety	7
2.1	Personnel and target group	7
2.2	Responsibility of the operator	8
2.3	Working in a safety conscious manner	8
2.4	Operating conditions	9
2.5	Media	9
2.6	Use.....	10
2.7	Directives and standards	11
2.8	Customer service and repair.....	11
2.9	Disposal	12
3	Technical data.....	13
3.1	Technical data N022.18	13
3.2	Technical data N026.3.18	15
3.3	Technical data N026.1.2.18	17
4	Product description	20
5	Transport.....	23
	General	23
6	Setup and connection	25
6.1	Preparing for commissioning	26
6.2	Perform commissioning	26
7	Operation	31
7.1	Information on switching the pump on and off	31
8	Servicing	32
8.1	Servicing schedule.....	33
8.2	Cleaning.....	34
8.3	Replacing diaphragm and valves.....	34
9	Spare parts and accessories.....	48
9.1	Spare parts	48
9.2	Accessories.....	49
10	Troubleshooting	50
11	Returns.....	54

1 About this document

1.1 Using the operating instructions

The operating instructions are part of the pump.

- In the event of uncertainties with regard to the content of the operating instructions, please contact the manufacturer (contact data: see *www.knf.com*). Please have the type and serial number of the pump ready.
- Read the operating instructions before you commission the pump.
- Only pass on the full and unchanged operating instructions to any subsequent owner.
- Keep the operating instructions within reach at all times.

1.2 Exclusion of liability

The manufacturer assumes no liability for damages and malfunctions resulting from failure to observe the operating instructions.

The manufacturer assumes no liability for damages and malfunctions resulting from changes or modifications to the device and improper handling.

The manufacturer assumes no liability for damages and malfunctions resulting from impermissible spare parts and accessories.

1.3 Symbols and markings

Warning notice

A notice that warns you of danger is located here.

Possible consequences of a failure to observe the warning notice are specified here. The signal word, e.g., Warning, indicates the danger level.

→ Measures for avoiding the danger and its consequences are specified here.

Danger levels

Signal word	Meaning	Consequences if not observed
DANGER	warns of immediate danger	Death or serious injury or serious damage will result.
WARNING	warns of possible danger	Death, serious injury or serious damage is possible.
CAUTION	warns of a possibly dangerous situation	Minor injury or damage is possible.
NOTICE	Warns of possible damage	Damage is possible.

Tab.1: Danger levels

Other notices and symbols

→ An activity to be carried out is specified here (a step).

1. The first step of an activity to be carried out is specified here.

Other sequentially numbered steps follow.

i This symbol indicates important information.

Explanation of pictograms

Pictogram	Meaning
	General warning symbol
	Warning of hot surface
	Warning of electrical voltage
	Warning of poisonous substances
	Warning of hand injuries through crushing
	Observe the operating instructions
	General mandatory sign
	Unplug mains plug
	Use foot protection
	Use hand protection
	WEEE Symbol for separate tracking of electrical and electronic devices. The use of this symbol means that this product must be disposed of with normal household waste.
	Recycling

Tab.2: Explanation of pictograms

2 Safety

i Observe the safety notices in Chapters 6 *Setup and connection* [▶ 25] and 7 *Operation* [▶ 31].

2.1 Personnel and target group

Personnel Ensure that only personnel who have received the appropriate training and instruction carry out work on the pumps. This applies in particular to commissioning and maintenance work.
 Make sure that the personnel have read and understood the operating instructions, particularly the chapter on safety.

Target group

Target group	Definition
User	Laboratory worker
Specialized personnel	Specialized personnel are personnel who - have relevant professional training in the field covered in the particular section of text; - have current knowledge of the field covered in the particular section of text.

Tab.3: Target group

Who-does-what matrix

Lifecycle phase	User	Specialized personnel
Transport		X
Setup	X	X
Preparing for commissioning	X	X
Commissioning	X	X
Operation	X	X
Servicing		X
Troubleshooting		X
Disposal		X

Tab.4: Who-does-what matrix

2.2 Responsibility of the operator

The pumps are produced in accordance with the generally recognized rules of engineering, as well as the occupational health, safety and accident prevention regulations. Nevertheless, dangers can arise during their use that lead to injuries to the user or third parties or to damage to the pump or other property.

Make sure that no hazardous situation, physical harm or impairment of the pump can occur.

Operating
parameter

Only operate and install the pump under the operating parameters and operating conditions described in Chapters 2.4 *Operating conditions* [► 9] and 3 *Technical data* [► 13].

2.3 Working in a safety conscious manner

Observe the regulations on accident prevention and safety during all work on the pumps and during operation.

Avoid contact with the pump heads and housing parts because the pump heats up during operation.

Make sure that the pump has been disconnected from the mains and is de-energized while work is carried out on it.

When connecting the pump to the electrical mains, observe the corresponding safety rules.

Do not expose any body parts to the vacuum.

Ensure that no hazards arise from gas flowing when gas connections are open, from the effects of noise or from hot, corrosive, dangerous and environmentally hazardous gases.

Ensure that EMC-compliant installation of the pump is guaranteed at all times, and that no hazardous situations can arise therefrom.

Avoid the release of hazardous, toxic, explosive, corrosive, harmful or environmentally hazardous gases or vapors, e.g. by using suitable laboratory equipment with fume cupboard and ventilation control.

2.4 Operating conditions

Do not use the pumps if they are not in technically perfect working order; pumps must be operated in accordance with their intended use, with attention to safety and hazard considerations, and in conformance with the operating instructions at all times.

The pumps must only be operated in the fully assembled and supplied condition.

Make sure that the installation location is dry and that the pump is protected from rain, splash water, gushing water, dripping water and other contamination.

Check the tightness of the connections between the pipes of the application and the pump (or the connection of the pump) at regular intervals. Leaky connections carry the risk of releasing dangerous gases and vapors from the pump system.

The components that are to be connected to the pumps must be designed according to the pneumatic data of the pumps.

2.5 Media

Requirements of pumped media	<p>Before transferring a medium, check whether the medium can be transferred danger-free in the specific application.</p> <p>Take note of any change in the state of matter (condensation, crystallization).</p> <p>Before using a medium, check the compatibility of the media-contacting components (see 3 <i>Technical data</i> [▶ 13]) with the medium.</p> <p>Only transfer gases that remain stable under the pressures and temperatures that arise in the pump.</p>
Handling of hazardous media	<p>Upon breakage of the diaphragm and/or leaks, the transferred medium mixes with the air in the surroundings and/or in the pump housing. Make sure that a dangerous situation cannot arise as a result.</p> <p>When pumping hazardous media, follow the safety regulations that apply for working with these media.</p>

- Handling of combustible media
- Note that the pump is not designed to be explosion-proof.
- Make certain that the temperature of the medium is always sufficiently below the ignition temperature of the medium so as to prevent ignition or explosion. This also applies for abnormal operating situations.
- At the same time, note that the temperature of the medium rises as the pump compresses the medium.
- Therefore, make certain that the temperature of the medium also remains sufficiently below the ignition temperature of the medium even when it is compressed to the maximum permissible operating pressure of the vacuum system. The maximum permissible operating pressure of the vacuum system is given in Chapter 3 *Technical data* [▶ 13].
- Make certain that the permissible ambient temperature (see 3 *Technical data* [▶ 13]) is not exceeded.
- Where applicable, also take into account external energy sources (such as radiated heat sources) that might heat the medium further.
- In case of doubt, contact KNF Customer Service.

2.6 Use

2.6.1 Proper use

The pumps are intended exclusively for delivering gases and vapors.

2.6.2 Foreseeable misuse

The pumps must not be operated in explosive atmospheres.

The pumps are not suitable for transferring the following:

- Dusts
- Liquids
- Aerosols
- Biological and microbiological substances
- Fuels
- Fibers

- Oxidizing agents
- Foodstuffs.

As standard, the pumps must not be used for simultaneous generation of a vacuum and positive pressure.

No overpressure may be applied to the suction side of the pump.

The pump must not be used if reactive explosive, or otherwise dangerous mixtures can occur (e.g. with the medium) when the gas ballast valve of the pump is open.

2.7 Directives and standards

EU/EC
Directives /
Standards

The pumps conform the following directives/ordinances:

- 2011/65/EU (RoHS)
- 2014/30/EU (EMC)
- 2006/42/EC (MD)

The following harmonized/listed standards are met:

- EN 61326-1 – Class A
- EN 1012-2
- EN ISO 12100
- EN 61010-1
- EN IEC 63000

Per IEC 664, the pumps comply with:

- Overvoltage category II
- Degree of soiling 2

2.8 Customer service and repair

Customer service and repairs

The pumps are maintenance-free. However, KNF recommends that the pumps be checked regularly for noticeable changes in noise and vibration.

Have repairs to the pumps carried out solely by the KNF customer service personnel responsible for such tasks.

Housings with electrically live components may only be opened by specialist personnel.

Use only genuine spare parts from KNF when performing servicing work.

2.9 Disposal

Environmental
protection
WEEE

Store the pump and all accessories in accordance with the environmental provisions. Observe the national and international regulations. This applies in particular to parts that are contaminated with toxic substances.

If you no longer need your packaging materials (e.g. for return shipment or other transport of the vacuum system), dispose of them in an environmentally friendly manner.

This product is marked in conformance with the EU directive on the disposal of waste electrical and electronic equipment (WEEE). Old devices must not be disposed of with household waste. Proper disposal and recycling help to protect natural resources and the environment. The end user is responsible for disposing of old devices according to the national and international regulations. Alternatively, KNF products (old devices) may also be returned to KNF for a fee (see chapter Returns).

3 Technical data

3.1 Technical data N022.18

Materials of media-contacting components

Assembly	Material AN	Material AT
Pump head	Aluminum	Aluminum
Diaphragm	CR	PTFE-coated
Valve	Stainless steel	Stainless steel
Seal	CR	FPM

Tab.5:

Pneumatic parameters

Parameter	Value N022AN.18	Value N022AT.18
Max. permissible operating pressure [bar rel]	4.0	4.0
Ultimate vacuum [mbar abs.]	100	100
Flow rate at atm. pressure [l/min]*	15.0	13.0

Tab.6: *Liters in the standard state based on ISO 8778 and ISO 21360-1/2 (1013 hPa, 20°C)

Electrical performance

Parameter	Value N022.18	Value N022.18	Value N022.18	Value N022.18
Voltage [V]	100	115	230	230
Frequency [Hz]	50/60	60	50	50
Power P ₁ [W]	140	130	100	120
Max. current draw [A]	2.0	1.6	0.7	1.0
Motor protection class (DIN EN 60529 / IEC 60529)	IP20	IP20	IP20	IP44
Max. permissible supply voltage fluctuations	± 10%	± 10%	± 10%	± 10%

Tab.7: Electrical performance

Thermal switch

i The motors are equipped standard with a thermal switch for monitoring against overheating.

Pneumatic connections

Pneumatic connections	Value
Hose connection [mm]	ID6

Tab.8: *Acc. to ISO 228

Other parameters

Parameter	Value
Permissible ambient temperature [°C]	+ 5 to + 40
Permissible media temperature [°C]	+ 5 to + 40
Relative air humidity	80% for temperatures to 31 °C, decreasing linearly to 50% at 40 °C (non-condensing).
Maximum installation altitude [m above sea level]	2000
Dimensions [mm] N022.18 IP20 N022.18 IP44	203 x 194 x 145 260 x 193 x 180

Tab.9:

Weight

Parameter	Value N022.18 IP20	Value N022.18 IP44
Weight [kg]	4.0	5.5

3.2 Technical data N026.3.18**Materials of media-contacting components**

Assembly	Material AN	Material AT
Pump head	Aluminum	Aluminum
Diaphragm	CR	PTFE-coated
Valve	Stainless steel	Stainless steel
Seal	CR	FPM

Tab.10:

Pneumatic parameters

Parameter	Value N026.3AN.18	Value N026.3AT.18
Max. permissible operating pressure [bar rel]	-	-
Ultimate vacuum [mbar abs.]	20	25
Flow rate at atm. pressure [l/min]*	22.0	18.0

Tab. 11: *Liters in the standard state based on ISO 8778 and ISO 21360-1/2 (1013 hPa, 20°C)

Electrical performance

Parameter	Value N026.3.18	Value N026.3.18	Value N026.3.18
Voltage [V]	100	115	230
Frequency [Hz]	50/60	60	50
Power P ₁ [W]	190	180	170
Max. current draw [A]	2.6	2.0	0.85
Motor protection class (DIN EN 60529 / IEC 60529)	IP20	IP20	IP20
Max. permissible supply voltage fluctuations	± 10%	± 10%	± 10%

Tab. 12: Electrical performance

Thermal switch

i The motors are equipped standard with a thermal switch for monitoring against overheating.

Pneumatic connections

Pneumatic connections	Value
Hose connection [mm]	ID9

Tab. 13: *Acc. to ISO 228

Other parameters

Parameter	Value
Permissible ambient temperature [°C]	+ 5 to + 40
Permissible media temperature [°C]	+ 5 to + 40
Relative air humidity	80% for temperatures to 31 °C, decreasing linearly to 50% at 40 °C (non-condensing).
Maximum installation altitude [m above sea level]	2000
Dimensions [mm]	243 x 192 x 185

Tab. 14:

Weight

Parameter	Value N026.3.18
Weight [kg]	6.3

3.3 Technical data N026.1.2.18**Materials of media-contacting components**

Assembly	Material AN	Material AT
Pump head	Aluminum	Aluminum
Diaphragm	CR	PTFE-coated
Valve	Stainless steel	Stainless steel
Seal	CR	FPM

Tab. 15:

Pneumatic parameters

Parameter	Value N026.1.2AN.18	Value N026.1.2AT.18
Max. permissible operating pressure [bar rel.]	2.0	2.0
Ultimate vacuum [mbar abs.]	100	100
Flow rate at atm. pressure [l/min]*	39.0	31.0

Tab. 16: *Liters in the standard state based on ISO 8778 and ISO 21360-1/2 (1013 hPa, 20°C)

Electrical performance

Parameter	Value N026.1.2.18	Value N026.1.2.18
Voltage [V]	115	230
Frequency [Hz]	60	50
Power P ₁ [W]	180	170
Max. current draw [A]	2.0	0.85
Motor protection class (DIN EN 60529 / IEC 60529)	IP20	IP20
Max. permissible supply voltage fluctuations	± 10%	± 10%

Tab. 17: Electrical performance

Thermal switch

i The motors are equipped standard with a thermal switch for monitoring against overheating.

Pneumatic connections

Pneumatic connections	Value
Hose connection [mm]	ID9

Tab. 18: *Acc. to ISO 228

Other parameters

Parameter	Value
Permissible ambient temperature [°C]	+ 5 to + 40
Permissible media temperature [°C]	+ 5 to + 40
Relative air humidity	80% for temperatures to 31 °C, decreasing linearly to 50% at 40 °C (non-condensing).
Maximum installation altitude [m above sea level]	2000
Dimensions [mm]	254 x 192 x 185

Tab. 19:

Weight

Parameter	Value
	N026.1.2.18
Weight [kg]	6.3

4 Product description

- 1 Outlet
- 2 Inlet
- 3 Power switch
- 4 Carrying handle


Fig.1: Design N022.18

- 1 Outlet
- 2 Inlet
- 3 Power switch
- 4 Pneumatic connection 1
- 5 Carrying handle
- 6 Pneumatic connection 2


Fig.2: Design N026.1.2.18

- 1 Outlet
- 2 Inlet
- 3 Power switch
- 4 Carrying handle
- 5 Pneumatic connection
- 6 Silencer / inlet filter (accessories)


Fig.3: Design N026.3.18

Function of a diaphragm pump

- 1 Outlet valve
- 2 Inlet valve
- 3 Transfer chamber
- 4 Diaphragm
- 5 Eccentric
- 6 Connecting rod


Fig.4: Function of a diaphragm pump

Diaphragm pumps transfer, compress (depending on the version) and evacuate gases and vapors.

The elastic diaphragm (4) is moved up and down by the eccentric (5) and the connecting rod (6). In the downwards stroke, it aspirates the gas to be transferred via the inlet valve (2). In the upwards stroke, the diaphragm presses the medium out of the pump head via the outlet valve (1). The transfer chamber (3) is separated from the pump drive by the diaphragm.

5 Transport

General

CAUTION

Personal injury and/or property damage due to incorrect or improper transport of the pump

In the event of incorrect or improper transport, the pump can fall down, be damaged or injure persons.

- Use suitable auxiliary means if necessary (carrying strap, lifting gear, etc.).
- Where appropriate, wear suitable personal protective equipment (e.g., safety shoes, safety gloves).

CAUTION

Risk of injury from sharp edges on the packaging

There is a risk of injury from cutting on the sharp edges when grabbing corners or when opening the packaging.

- Where appropriate, wear suitable personal protective equipment (e.g., safety shoes, safety gloves).

Only for two-headed pumps:

CAUTION

Personal injury and/or property damage due to incorrect or improper transport of the pump

If the pump is raised and/or carried on the connection during transport, leaks and/or damages to the pump may result.

→ Do not carry the pump on the pneumatic connection.

- Transport the pump in the original packaging to the installation location.
- Keep the original packaging of the pump (e.g. for later storage).
- Inspect the pump for transport damage after receiving it.
- Document any transport damage in writing.
- Remove any transport safeguards on the pump prior to commissioning.

Parameter

Parameter	Value
Storage temperature [°C]	+ 5 to + 40
Transport temperature [°C]	- 10 to + 60
Permissible humidity (non-condensing) [%]	30 to 85

Tab.20: Transport parameters and storage parameters

NOTICE

Prior to commissioning, make sure that the pump has reached the ambient temperature (3 *Technical data* [▶ 13]).

6 Setup and connection

- Only connect the pump in accordance with the operating parameters and conditions described in Chapter 3 *Technical data* [13].
- Observe the safety instructions (see Chapter Safety).
- Before connecting, store the pump at the installation location to allow it to reach the room temperature (no condensate may form).

Cooling air supply

WARNING

Danger of burning on hot surfaces
Hot surfaces could occur if the pump overheats.

- When installing the pump, make sure that sufficient cooling air in-feed and discharge is ensured.

Installation location

- Make sure that the installation location is dry and that the pump is protected against rain, splash, gushing, and drip water as well as from other contamination.
- Select a secure location (flat surface) for the pump.
- Protect the pump from dust.
- Protect the pump from vibration, impact and external damage.
- Make sure that it is easy to operate the power switch.

6.1 Preparing for commissioning

Before turning on the pump, make sure of the following points:

	Necessary operating requirements
Pump	- Connect all hoses correctly
Pump	- Data of the voltage supply system are consistent with the details on the type plate of the power supply. - Pump outlet not closed or restricted.
Pump	- Fan opening not blocked.

Tab.21: Operating requirements for commissioning

6.2 Perform commissioning

WARNING

Risk of burns from hot pump parts and/or hot medium

Some pump parts may be hot during or after operation of the pump.

- Allow the pump to cool after operation.
- Take protective measures to protect against touching hot parts.

Risk of injury from bursting hoses during pressure applications due to excessively high temperatures

When operating the pump in pressure applications, hoses that are not designed for the head temperatures of the pump at the respective operating point could become porous and burst.

- Use temperature-resistant pressure hoses at the pneumatic connections.
- Wear protective equipment if necessary (e.g., safety gloves, hearing protection).

Injury to eyes

Coming too close to the inlet/outlet of the pump may result in injury to the eyes due to the present vacuum/operating pressure.

- Do not look into the pump inlet/outlet during operation.

- Only operate the pump in accordance with the operating parameters and operating conditions described in Chapter 3 *Technical data* [▶ 13].
- Ensure the proper use of the pump (See Chapter 2.6.1 *Proper use* [▶ 10]).
- Eliminate the possibility of improper use of the pump (see Chapter *Improper use*).
- Note and follow the safety instructions (see chapter 2 *Safety* [▶ 7]).

Risk of pump head bursting due to excessive pressure increase

- Do not exceed the maximum permissible operating pressure (see 3 *Technical data* [▶ 13]).
- Monitor the pressure during operation.
- If the pressure exceeds the maximum permissible operating pressure of the pump: immediately switch off the pump and remedy the fault (see Chapter *Troubleshooting*).
- Only throttle or regulate the air or gas quantity on the suction line to prevent the maximum permissible operating pressure from being exceeded.
- If the air quantity or gas quantity on the pressure line is throttled or regulated, make sure that the maximum permissible operating pressure at the pump is not exceeded.
- Ensure that the pump outlet is not closed or restricted.

DANGER

Risk of dangerous gas mixtures during pump operation

Depending on the medium being transferred, breakage of the media-contacting components can result in a dangerous mixture if the medium mixes with the air in the compressor housing or the surroundings.

- Before using a medium, check the compatibility of the media-contacting components (see 3 *Technical data* [▶ 13]) with the medium.

NOTICE

Excessive pressure, with all of the associated hazards, can be prevented by means of a bypass line with a pressure relief valve between the pressure side and suction side of the pump. Further information is available from KNF Customer Service (contact data: see www.knf.com).

- Pump standstill → Establish normal atmospheric pressure in the lines while the pump is at a standstill (relieve pump pneumatically).

Connecting the pump

i A marking on the pump head indicates the flow direction.

1. Remove the protective caps from the pneumatic connections of the pump.
2. Mount suction filter or silencer accessory parts (if present).

**NOTICE**

If the pump is used as a vacuum pump, install a silencer on the pressure side if necessary.

If the pump is used as a compressor (not for pumps with .3 connection), install an inlet filter on the suction side if necessary.

-
- | | |
|----------------------|---|
| Connected components | 3. Connect the lines to the pneumatic inlet and outlet.
4. Only connect components to the pump that are designed for the pneumatic data of the pump (see Chapter 3 <i>Technical data</i> [▶ 13]). |
| Pump discharge | 5. When using as a vacuum pump: Safely drain the pump discharge at the pneumatic outlet of the pump.
6. Lay the line at the pneumatic inlet and the line at the pneumatic outlet in a downward slope so that no condensate can run into the pump.
7. Plug the plug of the power cable into a properly installed, grounded socket. |

7 Operation

7.1 Information on switching the pump on and off

Switching on the pump

NOTICE

The pump must not be started up against pressure or vacuum during switch-on. This also applies after a brief power interruption during operation. If a pump starts up against pressure or vacuum, the pump may block, causing the thermal switch to respond and the pump to switch off.

-
- Ensure that no pressure or vacuum is present in the lines when switching on.
 - Switch the pump on with the power switch.

Switching off/decommissioning the pump

- When transferring aggressive media, flush the pump before switching off to extend the service life of the diaphragm (see Chapter 8 *Servicing* [▶ 32]).
- Switch the pump off with the power switch.
- Establish normal atmospheric pressure in the lines (relieve pump pneumatically).
- Pull out power plug of the pump.

8 Servicing

NOTICE

Servicing the pump

Damage to the pumps can result from failure to observe the applicable legal regulations and procedures for the location or intervention by untrained or uninstructed personnel.

- Servicing may only be performed according to the legal regulations (e.g. work safety, environmental protection) and provisions.
 - Servicing may only be performed by specialized personnel or trained and instructed personnel.
-

8.1 Servicing schedule

Risk of injury when not using original parts

Failure to use original parts will result in a loss of pump functionality and safety.

The validity of the CE conformity is rendered void if genuine parts are not used.

- Use only genuine spare parts from KNF when performing maintenance/repair work.

Component	Servicing interval
Pump	<ul style="list-style-type: none"> → Perform periodic inspections for external damage or leakage. → Periodically check for noticeable changes to noises and vibrations.
Diaphragms and valves	<ul style="list-style-type: none"> → At the latest, replace when the performance decreases.
Inlet filter (accessories)	<ul style="list-style-type: none"> → Replace if soiled
Silencer (accessories)	<ul style="list-style-type: none"> → Replace if soiled

Tab.22: Servicing schedule

8.2 Cleaning

NOTICE

During cleaning work, ensure that no fluids enter the interior of the housing.

8.2.1 Flushing the pump

- Before switching off, flush the pump with air at atmospheric conditions (ambient pressure) for about 5 minutes (if necessary for safety reasons: with an inert gas).

8.2.2 Cleaning the pump

- Only clean the pump with a damp cloth and non-flammable cleaning agents.
- If compressed air is present, blow out the parts.

8.3 Replacing diaphragm and valves

8.3.1 Change diaphragm and valves (A_ version)

Spare parts

Spare part*	Item number**	Number per pump head
Diaphragm	(F)	1
Countersunk screw	(D)	1
Valve spring	(M, P)	2
Seal	(V)	1

Tab.23: Spare parts

*According to spare parts list, Chapter 9.1 Spare parts [▶ 48]

** According to Fig. 5

Tool	Quantity	Tool/material
	1	Allen key, 3 mm
	1	Allen key, 4 mm
	1	Screwdriver blade width 6.5
	1	Screwdriver blade width 4.0
	1	Socket wrench 5.5 mm
	1	Pencil

Tab.24: Tool

Information on the procedure The diaphragm and valves are the only wearing parts in the pumps. They are easy to replace.

The valves and diaphragm should generally be replaced at the same time. If the diaphragm is not replaced at the same time as the valves, the specified output of the pump can no longer be ensured after the maintenance is performed.

With multi-headed pumps, parts may be interchanged between the individual pump heads.

→ Replace the diaphragm, valve springs and seals of the individual pump heads in sequence.

WARNING

Health hazard due to dangerous substances in the pump

Depending on the medium being transferred, caustic burns or poisoning is possible.

→ Wear protective equipment if necessary, e.g., protective gloves, goggles.

→ Clean the pump with suitable measures.

CAUTION

Risk of burns from hot pump parts

The pump head or motor may still be hot after operation of the pump.

→ Allow the pump to cool after operation.

- A** Housing
- B** Allen screw
- C** Membrane head
- D** Countersunk screw
- E** Retainer plate
- F** Diaphragm
- G** Allen screw
- H** Cover
- I** Flywheel
- K** Connecting rod
- M** Valve spring
- P** Valve spring
- S** Allen screw
- T** Cover
- U** Nut
- V** Seal
- W** Screw
- X** Washer


Fig.5: Pump parts

1. For pumps N026.1.2.18:
 Pull off the hose on a pump head at the suction-side pneumatic connection between the pump heads.
 Loosen the hose clamp on the pressure-side pneumatic connection on a pump head and pull off the hose.

2. For pumps N026.3.18:
Pull off the hose of the pneumatic connection at a pump head.
3. Mark the position of the diaphragm head (C) with respect to the housing (A) with a pencil mark.
4. Loosen the 4 Allen screws (B) and remove the diaphragm head (C).
5. Loosen the countersunk screw (D).
6. Remove the retainer plate (E) and the diaphragm (F).
7. Loosen the 4 cheese head screws (G) and remove the cover (H).
8. Turn the flywheel (I) until connecting rod (K) is in the center position.
9. Put on the new diaphragm (F).
10. Place the retainer plate (E) on the diaphragm (F) and tighten it with the new countersunk screw (D) (tightening torque: 5.0 Nm).

i The self-locking countersunk screw (D) can only be used once.

11. Change the valve springs (M) and (P):
 - Loosen the Allen screws (S).
 - Lift off the cover (T) and the seal (V).
 - Loosen the nut (U) with a socket wrench.
 - Pull out the screw (W).
 - Remove the valve springs (P) and (M).
 - Fix the new valve springs (P) and (M) with the screw (W), washer (X) and nut (U).
 - Put on the cover (T) with new seal (V).
 - Tighten the Allen screws (S).
12. Place the diaphragm head (C) in the installation position marked with a pencil and tighten the Allen screws (B) evenly crosswise (tightening torque: AN and AV version: 6.5 Nm; AT version: 5.5 Nm)
13. Check for smooth operation when turning the flywheel (J).

14. For two-headed pumps:
Perform work steps 3 to 13 on the second pump head.
15. Fasten the cover (H) with the 4 cheese head screws (G).
16. For two-headed pumps:
Pull the hose (pumps N026.1.2.18: the hoses) of the pneumatic connection back onto the hose connector.
For pump types N026.1.2.18: Retighten the hose clamp on the pneumatic connection on the pressure side.

8.3.2 Change diaphragm and valves (ST version on project basis)

Spare parts

Spare part*	Item number**	Number per pump head
Diaphragm	(F)	1
Countersunk screw	(D)	1
Valve plate	(Z)	1

Tab.25: Spare parts

*According to spare parts list, Chapter 9.1 Spare parts [▶ 48]

** According to Fig. 6

Tool

Quantity	Tool/material
1	Allen key, 4 mm
1	Screw driver blade width 6.5
1	Adjustable face spanner wrench for two-hole nuts or KNF wrench for retainer plate (.9 version only)
1	Pencil

Tab.26: Tool

Information on the procedure

The diaphragm and valves are the only wearing parts in the pumps. They are easy to replace.

The valves and diaphragm should generally be replaced at the same time. If the diaphragm is not replaced at the same time as the valves, the specified output of the pump can no longer be ensured after the maintenance is performed.

With multi-headed pumps, parts may be interchanged between the individual pump heads.

→ Replace the diaphragm, valve springs and seals of the individual pump heads in sequence.

**WARNING**

Health hazard due to dangerous substances in the pump

Depending on the medium being transferred, caustic burns or poisoning is possible.

- Wear protective equipment if necessary, e.g., protective gloves, goggles.
- Clean the pump with suitable measures.

**CAUTION**

Risk of burns from hot pump parts

The pump head or motor may still be hot after operation of the pump.

- Allow the pump to cool after operation.

- A** Housing
- D** Countersunk screw
- E** Retainer plate
- F** Diaphragm
- G** Allen screw
- H** Cover
- I** Flywheel
- K** Connecting rod
- W** Membrane head
- X** Intermediate plate
- Y** Allen screw
- Z** Valve plate


Fig.6: Pump parts

17. For pumps N026.1.18 and N026.3.18:
Pull off the hose of the pneumatic connection at a pump head.
18. For pumps N026.2.18:
Loosen the hose clamp of the pneumatic connection on a pump head and pull off the hose.

19. Mark the position of the diaphragm head (W) and the intermediate plate (X) with respect to the housing (A) with a pencil mark.
 20. For all pumps except .9 version:
Loosen the countersunk screw (D) and remove the retainer plate (E) and the diaphragm (F).
 21. For .9 version:
Use wrench for retainer plate to loosen the retainer plate (E) counterclockwise; remove the retainer plate (E) and the diaphragm (F).
 22. Loosen the 4 cheese head screws (G) and remove the cover (H).
 23. Turn the flywheel (I) until connecting rod (K) is in the center position.
 24. Put on the new diaphragm (F).
 25. For all pumps except .9 version:
Place the retainer plate (E) on the diaphragm (F) and tighten it with the new countersunk screw (D) (tightening torque: 5.0 Nm).
- i** The self-locking countersunk screw (D) can only be used once.
26. For .9 version:
Place the retainer plate (E) on the diaphragm (F) and tighten it clockwise with the wrench for the retainer plate (tightening torque: 5.0 Nm).
 27. Place the intermediate plate (X) on the diaphragm (F) according to the pencil marking on the housing.
 28. Place the new valve plate (Z) on the intermediate plate (X) (for alignment, see Fig. 6).
 29. Place the diaphragm head (W) in the installation position marked with a pencil and tighten the Allen screws (Y) evenly crosswise (tightening torque: 5.5 Nm)
 30. Check for smooth operation when turning the flywheel (I).
 31. For two-headed pumps:
Perform work steps 3 to 14 on the second pump head.
 32. Fasten the cover (H) with the 4 cheese head screws (G).

33. For two-headed pumps:
 Pull the hose of the pneumatic head connection back onto the hose connector.
 For pump types N026.2.18: Retighten the hose clamp on the pneumatic connection on the pressure side.

8.3.3 Change diaphragm and valves (SP version on project basis)

- Requirements → Disconnect the motor from the mains and ensure that it is voltage-free.
 → Free the pump of hazardous materials.

Spare part*	Item designation**	Number per pump head
Diaphragm	(F)	1
Valve spring (stainless steel)	(M)	2
O-ring	(B)	2

Tab.27: *According to spare parts list, Chapter Spare parts

** According to Fig. 7

Tools and material

Quantity	Tool/material
1	Allen key, 4 mm
1	Screw driver blade width 4.5 mm
1	Adjustable face spanner wrench for nuts with two holes or KNF wrench for retainer plate (see Accessories).
1	Pencil

Tab.28:

- Information on the procedure → Always replace diaphragm, valve springs and O-rings together to maintain the performance of the pump.
 With multi-headed pumps, parts may be interchanged between the individual pump heads.
 → Replace the diaphragm and valve springs of the individual pump heads in sequence.

**WARNING**

Health hazard due to dangerous substances in the pump

Depending on the medium being transferred, caustic burns or poisoning is possible.

- Wear protective equipment if necessary, e.g., protective gloves, goggles.
- Clean the pump with suitable measures.

Work steps

i The following item numbers refer to Fig. 7.

1. Only for two-headed pumps:
Open the pneumatic connection between the two heads; loosen one of the union nuts to do this; be careful not to kink the connecting hose.
2. Removing the housing cover (H):
Loosen the fastening screws (G) of the housing cover (H) and remove the housing cover.
3. On a pump head, mark housing (A), intermediate plate (D) and head plate (C) with a continuous line made with a pencil. This helps to avoid incorrect assembly later.
4. Loosen the four hexagon socket head cap screws (Y) and remove the head plate (C) and the intermediate plate (D) from the housing.
5. Remove the O-rings (B) from the head plate (C).
6. Loosen the valve fastening screws (W) in the head plate and the intermediate plate and remove the valve springs (M).
7. For all pumps except .9 version:
Loosen the countersunk screw (L) and remove the retainer plate (E) and the diaphragm (F).

8. For .9 version:
Use wrench for retainer plate to loosen the retainer plate (E) counterclockwise; remove the retainer plate (E) and the diaphragm (F).
 9. Check all parts for contamination and clean the parts if necessary.
 10. Mounting the new valve springs (M):
 11. Check the pretension of the valve springs; to do this, hold the valve springs lengthwise between two fingers and press the valve springs together very lightly.
 12. Place the valve springs (M) on the valve seat of the head plate (C). The slight curvature of the valve spring caused by the pretension must face the head plate (C).
 13. Tighten the valve fastening screw (W) with a screw driver.
 14. Place the valve spring (M) on the valve seat of the intermediate plate (D). The slight curvature of the valve spring caused by the pretension must face the intermediate plate (D).
 15. Place the washer (X) on the valve fastening screw (W). Then tighten the valve fastening screw (W) with a screw driver.
 16. Check whether the valve springs are properly seated.
 17. Insert the new O-rings (B) in the head plate (C).
 18. Mount the diaphragm with the retainer plate.
 19. For all pumps except .9 version:
Place the retainer plate (E) on the diaphragm (F) and tighten it with the new countersunk screw (L) (tightening torque: 5.0 Nm).
- i** The self-locking countersunk screw (L) can only be used once.
20. For .9 version:
Place the retainer plate (E) on the diaphragm (F) and tighten it clockwise with the wrench for the retainer plate (tightening torque: 5.0 Nm).
 21. Turn the flywheel (I) to move the connecting rod (K) to the middle position.

22. Place the new diaphragm with the retainer plate on the connecting rod (K).
23. Screw in the retainer plate (E) clockwise with the wrench for retainer plate and hand tighten.
24. Place the intermediate plate (D) on the housing according to the pencil marking.
25. Place the head plate (C) on the intermediate plate (D) according to the pencil marking.
26. Uniformly tighten the hexagon socket head cap screws (Y) crosswise.
27. Check the pump for smooth running by turning the flywheel (I).
28. For two-headed pumps:
Perform work steps 3 to 25 on the second head.
29. Remount the housing cover (H).
30. For two-headed pumps only: Remount the pneumatic connection.
31. Checking the pump head (pump heads) and the pneumatic connections for leaks:

**NOTICE**

To ensure the required gas tightness of the pump following servicing, a leak test is to be performed.

**WARNING**

Risk of injury and poisoning from leaks

- Before recommissioning the pump, check the pump heads and pneumatic connections for leaks. Leaks may lead to poisoning, chemical burns or similar injuries.
-

- A** Housing
- B** O-ring
- C** Head plate
- D** Intermediate plate
- E** Retainer plate
- F** Diaphragm
- G** Fastening screw
- H** Housing cover
- I** Flywheel
- K** Connecting rod
- L** Countersunk screw
- M** Valve spring (stainless steel)
- W** Valve fastening screw
- X** Washer
- Y** Hexagon socket head cap screw


Fig.7: Pump parts for pump models with smooth heads and valve springs (stainless steel, symbolic)

9 Spare parts and accessories

i To order spare parts and accessories, please contact your KNF sales partner or KNF Customer Service (contact data: see www.knf.com).

9.1 Spare parts

Spare part set

Spare part set	Order number
N022AN.18	030305
N022AT.18	032490

Tab.29: Spare part set

Spare part set	Order number
N026.1.2AN.18	118951
N026.1.2AT.18	118953

Tab.30: Spare part set

Spare part set	Order number
N026.3AN.18	118951
N026.3AT.18	118953

Tab.31: Spare part set

A spare part set consists of:

Spare parts	Item number*	Quantity
Diaphragm	(F)	2
Countersunk screw	(D)	2
Valve spring	(M, P)	4
Seals	(V)	2

Tab.32: Spare parts
Acc. to Figure_Pump parts

9.2 Accessories

N022.18

Accessories	For pump type	Order number
Silencer / inlet filter (G 1/4)	N022.18	000346
Overpressure relief valve 4 bar	N022AN.18	000351
Fine adjustment head with pressure gauge on pressure side	N022AN.18	000349
Fine adjustment head with vacuum gauge on suction side	N022AN.18	000350

Tab.33: Accessories

N026.18

Accessories	For pump type	Order number
Silencer / inlet filter (G 1/4)	N026.18	000352
Overpressure relief valve 2 bar	N026.1.2AN.18	003074
Fine adjustment head with pressure gauge on pressure side	N026.1.2AN.18	011867
Fine adjustment head with vacuum gauge on suction side	N026.1.2AN.18	011868

Tab.34: Accessories

10 Troubleshooting

Danger: electric shock can be life-threatening.

- All work on the pump may only be performed by an authorized specialist.
- Before working on the pump: Disconnect the pump from the power supply.
- Check and ensure that no voltage is present.

→ Check the pump (see following tables).

Pump not delivering	
Cause	Troubleshooting
No voltage in the electrical mains.	→ Check the circuit breaker for the room and switch it on if necessary.
Thermal switch or overtemperature protection of the pump has tripped	<ul style="list-style-type: none"> → Disconnect the pump from electrical mains. → Allow the pump to cool down. → Determine the cause of the overheating and rectify.
Connections or lines are blocked.	<ul style="list-style-type: none"> → Check the connections and lines. → Remove the blockage.
External valve is closed or filter clogged.	→ Check external valves and filters.
Condensation has collected in the pump head.	<ul style="list-style-type: none"> → Separate the source of the condensation from the pump. → Flush the pump with air at atmospheric pressure for a few minutes (if necessary for safety reasons: with an inert gas).
Diaphragm or valves are worn.	→ Replace the diaphragms and the valve plates/seals (see Chapter 8 Servicing [▶ 32]).

Tab.35: Troubleshooting: Pump not delivering

Flow rate, pressure or vacuum too low

The pump does not reach the performance stated in the technical data or data sheet.

Cause	Troubleshooting
Condensation has collected in the pump head.	<ul style="list-style-type: none"> → Separate the source of the condensation from the pump. → Flush the pump with air at atmospheric pressure for a few minutes (if necessary for safety reasons: with an inert gas). → If present, open the gas ballast and flush the pump head.
There is overpressure on the pressure side and at the same time vacuum or pressure above atmospheric pressure on the suction side.	<ul style="list-style-type: none"> → Change the pneumatic conditions.
Pneumatic lines or connection parts have insufficient cross-sections or are constricted.	<ul style="list-style-type: none"> → Disconnect the pump from the system to determine the output values. → Eliminate any constriction (e.g. valve). → Use lines or connection parts with a larger cross-section if necessary.
Leaks occur at connections, lines or pump head.	<ul style="list-style-type: none"> → Ensure the correct seating of the hoses on the hose connectors. → Ensure that the connections are correctly mounted. → Replace the leaky hoses. → Eliminate the leaks.
Connections or lines are completely or partially clogged.	<ul style="list-style-type: none"> → Check the connections and lines. → Remove any parts or particles that are causing blockages.
Head parts are soiled.	<ul style="list-style-type: none"> → Clean the head components.
Diaphragm or valves are worn.	<ul style="list-style-type: none"> → Replace the diaphragms and the valve plates/seals (see Chapter 8 <i>Servicing</i> [▶ 32]).

Flow rate, pressure or vacuum too low	
The pump does not reach the performance stated in the technical data or data sheet.	
Cause	Troubleshooting
Replaced diaphragm and valves.	<ul style="list-style-type: none"> → Check the pneumatic connection and hoses for leaks. → If necessary, carefully tighten the screws (B) and/or (S) (see Fig. 5) crosswise.

Tab.36: Troubleshooting: Flow rate, pressure or vacuum too low

Pump is switched on and not running; power switch is not illuminated	
Cause	Troubleshooting
Pump is not connected to the electrical mains.	→ Connect the pump to the electrical mains.
No voltage in the electrical mains.	→ Check the circuit breaker for the room and switch it on if necessary.

Tab.37: Troubleshooting: Pump is switched on and not running; power switch is not illuminated

Pump is switched on and not running; power switch is illuminated	
Cause	Troubleshooting
Pump has overheated, thermal switch has tripped.	<ul style="list-style-type: none"> → Pull power cable of the pump out of the socket. → Allow the pump to cool down. → Determine the cause of the overheating and rectify.

Tab.38: Troubleshooting: Pump is switched on and not running; power switch is illuminated

Fault cannot be rectified

If you are unable to identify any of the specified causes, send the pump to KNF Customer Service (contact data: see www.knf.com).

1. Flush the pump with air for a few minutes (if necessary for safety reasons: with inert gas) at atmospheric pressure to free the pump head of dangerous or aggressive gases (see Chapter Flushing the pump).
2. Clean the pump (see Chapter Cleaning the pump).
3. Send the pump together with completed Health and Safety Clearance and Decontamination Form to KNF, specifying the pumped medium.

11 Returns

Preparing for return

1. Flush the pump with air for a few minutes (if necessary for safety reasons: with inert gas) at atmospheric pressure to free the pump head of dangerous or aggressive gases (see Chapter 8.2.1 *Flushing the pump* [▶ 34]).
i Please contact your KNF sales partner if the pump cannot be flushed due to damage.
2. Remove the pump.
3. Clean the pump (see Chapter 8.2.2 *Cleaning the pump* [▶ 34]).
4. Send the pump together with the completed Health and Safety Clearance and Decontamination Form to KNF, stating the nature of the transferred medium.
5. Pack the device securely to prevent further damage to the product. If necessary, request original packaging for a fee.

Returns

KNF shall undertake to repair the pump only under the condition that the customer presents a certificate regarding the medium that is pumped and the cleaning of the pump. In this case too, old devices can be returned. Please follow the instructions at [knf.com/repairs here](http://knf.com/repairs).

Contact your KNF sales partner directly if you require additional support for your return service.

KNF Neuberger GmbH
Alter Weg 3
79112 Freiburg
Germany
Tel. +49 (0)7664/5909-0

E-mail: info.de@knf.com
www.knf.com

KNF worldwide

You can find our local KNF partners at: www.knf.com

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