



America

CERTIFICATE

No. U8 045055 0021 Rev. 00


Holder of Certificate: **KNF Neuberger GmbH**



Alter Weg 3
79112 Freiburg
GERMANY

Production Facility(ies): 045055

Certification Mark:



Product: **Laboratory Equipment**
(Vacuum pump for laboratory use)

Model(s): **Laboport N96**

Parameters:
Rated voltage: 24 VDC
Rated current: 0.9 A
Rated power: 19 W
Protection class: III
Rated max. pressure: +250 kPa (2.5 bar)

Tested according to: UL 61010-1:2012/R:2018-11
CAN/CSA-C22.2 No. 61010-1:2012/A1:2018-11

This product was voluntarily tested to the relevant safety requirements referenced on this certificate. It can be marked with the certification mark above. The mark must not be altered in any way. This product certification system operated by TÜV SÜD America Inc. most closely resembles system 3 as defined in ISO/IEC 17067. Certification is based on the TÜV SÜD "Testing and Certification Regulations". TÜV SÜD America Inc. is an OSHA recognized NRTL and a Standards Council of Canada accredited Certification body.

Test report no.: 028-713175319-000

Date, 2020-04-28

(Abdul Sabbagh)