

DOSING PUMP FOR LIQUIDS

SIMDOS 10



Hydraulic data

Flow rate	Dosing volume	Suction height	Max pressure
1 - 100 ml/min	1 - 999 ml	3 mWg	6 bar

Chemical resistance

Code	Pump head	Diaphragm	Valves	Filter ¹⁾
KT	PP	PTFE covered	FFKM Kalrez [®]	PEEK
TT	PVDF	PTFE covered	FFKM Kalrez [®]	PVDF
FT	PTFE	PTFE covered	FFKM Kalrez [®] *	PEEK
ST	STAINLESS STEEL	PTFE covered	FFKM Kalrez [®]	PEEK

1) *On request FFKM Chemraz® material for housing and filter mesh

Functions

S Version	RC Version		RCP Version
Manual operation	 Manual operation Analog control: 0-5 V, 0-10 V, 0-20 mA, 4-20 mA, from 1 to 100% Start/Stop through logic control (TTL) 	 Reset/Prime through logic control (TTL) Output signal fault alarm Cable for external control included 	Same as RCPlus RS 232

Technical data

- Main supply 100 240 V / 50 60 Hz
- Allowed ambient temperature +5 to +40° C
- Allowed liquid temperature +5 to +80° C
- Maximum viscosity allowed 150 cSt
- Accuracy +/-2% (nominal value)
- Repeatability +/-1%

- Power consumption 24 W
- Protection class IP65
- Dimensions: 150 x 93 x 144 mm
- Weight: 0.9 kg
- Hydraulic connection: compression fittings for 4/6 mm tube

ACCESSORIES

Valves Kit

	ID-Nr.
Chemraz [®] special valves for certain chemicals	168037





Diaphragm pressure control valve with support when the feeding tank is above the pump or to dose into vacuum

	ID-Nr.
Kit FDV 30 KTZ	166283
Kit FDV 30 TTZ	166284
Kit FDV 30 FTZ	166285



Filter

Тур	Material ¹⁾	Connector	ID-Nr.
FS 25 T	PVDF	for tubes ID 3.2 and 4 mm	165211
FS 25 X	PEEK	for tubes ID 3.2 and 4 mm	165213



Pulsation damper only use in continuous transfer mode from 30 to 100 ml/min

	ID-Nr.
FPD 10 KT SIMDOS	167817
FPD 10 TT SIMDOS	167818
FPD 10 FT SIMDOS	167819



Mounting plate

	ID-Nr.
for wall mounting	160473



Standfitting

	ID-Nr.
for retort stand	160474



Footswitch

	ID-Nr.
for impulse start and stop	155872



| KNF Flodos AG, reserves the right to mak



¹⁾ material for housing and filter mesh